

RAPID, FULLY AUTOMATED BIOMARKER TESTING IDYLLA™




Idylla™ A revolutionary,
fully automated system that makes
molecular testing convenient and fast.
Suitable for any lab.



IDYLLA™ – MAKING MOLECULAR TESTING CONVENIENT AND SUITABLE FOR ANY LAB

Idylla™, a fully automated, one step sample-to-result PCR based molecular diagnostics system that can be used in any laboratory setting thanks to its compact, scalable design and outstanding ease of use.



IDYLLA™ - A REVOLUTIONARY SYSTEM FOR PCR BIOMARKER TESTING IN ONCOLOGY AND INFECTIOUS DISEASE

IDYLLA™ TEST MENU



Oncology: diagnostic tests for solid and liquid biopsies in colorectal cancer, lung cancer and melanoma



Infectious diseases: diagnostic tests for COVID-19¹ and sepsis²



IDYLLA™ - A REVOLUTIONARY WORKFLOW



Obtain sample



Transfer sample to cartridge



Insert cartridge and start test

Just add your sample. Idylla™ will take care of the rest. Our fully automated system enables result reporting within 65 to 160 minutes depending on the assay.

(1) The Idylla™ SARS-CoV-2 Test is a CE-marked IVD in Europe validated for the qualitative detection of SARS-CoV-2 RNA. The Idylla™ SARS-CoV-2 Test has been validated by Biocartis and is filed for Emergency Use Authorization with the FDA. The FDA's independent review of this validation is pending. The Idylla™ SARS-CoV-2 Test uses SuperScript™ III. The SuperScript™ III trademark is owned by Life Technologies Corporation.

(2) SeptiCyt® RAPID is a CE-marked IVD, developed by Immunexpress Inc in collaboration with Biocartis. Biocartis has the exclusive distribution rights for the EU. The test is not available in all countries. Availability to be checked with local Biocartis representative.

AUTHORIZED
DISTRIBUTOR



MMJ BIOSYSTEMS PHILIPPINES, INC.

Telephone: (+632) 8808-2219; (+632) 8636-1485

Email: info@mmjbiosystems.com

www.mmjbiosystems.com

Biocartis NV
Generaal de Wittelaan 11B
2800 Mechelen - Belgium
+32 15 632 888
www.biocartis.com
customerservice@biocartis.com

Follow us on



The Idylla™ platform is a CE-marked IVD in Europe.

Idylla™ is available for sale in EU, USA and some other countries. Please check availability with the local Biocartis sales representative.

Biocartis and Idylla are registered trademarks in Europe, the United States and other countries.

The Biocartis and the Idylla trademarks and logos are used trademarks owned by Biocartis.

© Biocartis, October 2020

